


ICE Devices

Electrolysis Apparatus



This device allows one to easily set up an electrolysis demonstration on an overhead projector. While the actual chemistry might be too complicated for younger students, the phenomenon of electrolysis and conductivity of water could be shared with children of all ages. The evolution of gas bubbles and the different colors of an indicator produced at the anode and the cathode are a fascinating sight.

Disclaimer. This publication is intended for use by teachers or other adults or by students with proper supervision. Users are expected to follow normal safety precautions when performing the activities described in this manual. Neither the authors nor the publisher assume liability for the use of information in this publication.

Copyright © 1994, 2015 by University of Wisconsin Board of Regents; ICE Publication: *ICE Devices*

Creative Commons License Deed. Attribution-NonCommercial-ShareAlike 3.0 Unported [<http://creativecommons.org/licenses/by-nc-sa/3.0/>] A portion of this document is below; see complete license online.

You are free to

Share—copy and redistribute this material in any medium or format;

Adapt—remix, transform, and build upon the material.

You may do so as long as you follow the license terms:

- **Attribution**—You must give appropriate credit, provide a link to the license, and indicate if changes were made. You may do so in any reasonable manner, but not in any way that suggests the licensor endorses you or your use.
 - **Noncommercial**—You may not use the material for commercial purposes.
 - **Share Alike**—If you remix, transform, or build upon the material, you must distribute your contributions under the same license as the original.
- No additional restrictions**—You may not apply legal terms or technological measures that legally restrict others from doing anything the license permits.

Permission to copy: You have our permission to make as many printed copies as are needed for students in your classes of any or all of these activities.

ICE Institute for Chemical Education

Contact Information

ICE, Institute for Chemical Education

University of Wisconsin–Madison
Department of Chemistry
1101 University Avenue
Madison, WI 53706

toll free 888-220-9822 (toll-free number, U.S.)
telephone 608/262-3033
fax 608/265-8094
email ice@chem.wisc.edu and iceorders@chem.wisc.edu
(sorry, we cannot accept credit card orders by email)

Visit the ICE Website

Information about all ICE Kits and Publications as well as our Order Form may be found on our web site at <http://ice.chem.wisc.edu>. The site also has information about our programs in Outreach, Education, and Research as well as the Nanoscience theme that runs through what we do and what we publish.

Electrolysis Apparatus


By William Fausey

This device allows one to easily set up an electrolysis demonstration on an overhead projector. While the actual chemistry might be too complicated for younger students, the phenomenon of electrolysis and conductivity of water could be shared with children of all ages. The evolution of gas bubbles and the different colors of an indicator produced at the anode and cathode are a fascinating sight.

The electrolysis of a KI solution has been written up by Ken and Doris Kolb [*J. Chem. Educ.* 1986, 63, 517]. Immediately upon immersion of the electrodes, there is a copious flow of hydrogen gas bubbles at the cathode and within a few seconds a yellow–brown iodine coloration at the anode. After 30–60 seconds a drop of phenolphthalein placed near the cathode turns that area of the solution bright pink. A squirt of 1% starch solution near the anode gives that area the blue-black color of the starch-iodine complex.

Electrolysis of NaNO_3 using various indicators was described by J. Skinner [*J. Chem. Educ.* 1981, 58, 1017]. Any dilute aqueous solution, e.g., simple tap water, can be electrolyzed in this way. The anode and cathode can be identified by adding a few drops of red cabbage juice indicator near each electrode.

The instructions here call for construction from a wooden block. It has been suggested that this device can also be made from a styrofoam block, such as used by florists to hold flower arrangements. However, in that case the apparatus may not stand up on its own, but will have to be held in your hand.



A completed electrolysis apparatus.

Materials

- ___ 2 × 4 wood block (6 cm × 3.5 cm × 8.5 cm)
- ___ 2 cm × 0.5 cm × 8.5 cm wood strip, e.g., a tongue depressor or a paint stirrer cut to size
- ___ saw
- ___ hammer and small nail or awl
- ___ vise or C-clamps and work table
- ___ drill
- ___ 9/16"-dia., 1/2"-dia. and 1/4"-dia. drill bits
- ___ 1/4"-dia. 2"-long machine bolt
- ___ wing nut for 1/4"-dia. bolt
- ___ #5 or #6 self-starting screw, 1/2 to 3/8" long
- ___ 9V battery
- ___ 9V battery holder (e. g. Radio Shack #270-326 or from an electronics supplier)
- ___ 9V battery snap connector (e.g. Radio Shack #270-325 or from an electronics supplier)
- ___ 2 ea. alligator clips (Radio Shack #270-375 or from an electronics supplier)
- ___ pliers
- ___ 2 ea. sharpened pencils
- ___ petri dish [e. g. Flinn Scientific AP8170; also widely available online from lab supply houses]
- ___ red cabbage juice (keep refrigerated)
- ___ 1 cup distilled water
- ___ blender
- ___ sieve
- ___ 2-cup mixing bowl or beaker

Procedure

Preparing the electrode holder

1. Drill a 1/4" hole in the center of the wood strip (Figure 1). This strip will serve as the electrode (pencil) holder.




Figure 1. Drilling the wood strip.

- Using a nail or awl start a small hole in the center of the $3.5\text{ cm} \times 8.5\text{ cm}$ surface of the wood block (Figure 2).




Figure 2. Making the hole for the battery holder.

- Cut a wedge from the top of the wood block as shown in Figure 3.
- Drill a $1/2''$ hole perpendicular to the sloped surface all the way through the block. Countersink the hole (that is, drill the hole out from the other side) about half way using a $9/16''$ drill bit.




Figure 3. Cutting and drilling the wood block.

- Attach the battery holder to the wood block, fixing it in place with the #5 or #6 self-starting screw (Figure 2).

- Near the eraser end of each of the pencils, remove enough wood to expose the lead.
- Insert the 1/4" machine bolt from the bottom of the block. Place the pencils on the slanted surface with their sharpened ends pointing down. Insert the machine bolt into the hole of the wood strip, thread on the wing nut and tighten it down so that the pencils do not budge. The pencils should be fairly parallel to each other and about 3 cm apart.
- Using the pliers, crimp the alligator clips to the wires on the battery clip.
- Clamp the battery in the battery holder and attach the battery clip. Attach the alligator clips to the exposed pencil leads.

Making Cabbage Juice


- Coarsely chop the red cabbage.
- Place the cabbage into a blender or food processor. Add about 250 mL (1 cup) of distilled water and blend the mixture until the cabbage has been pulverized.
- Pour the mixture through a sieve held over a container to catch the liquid. This strained liquid is the cabbage juice indicator. (**Note:** If left unrefrigerated, cabbage juice will spoil and become *very* smelly.)

Presentation

- Place the petri dish with the solution that you will be electrolyzing on the overhead projector.
- Place the electrode apparatus next to the petri dish so that the pencil leads are in the solution. Observe the electrolysis.
- Place several drops of the red cabbage juice indicator into the solution near the electrodes.
- Observe what happens near each electrode. Determine which electrode is the cathode and which one is the anode.

The actual reactions at the electrodes are:

Electrolysis of KI(aq)



Electrolysis of NaNO₃(aq)

